Module Interface Specification for CXR

Team 27, Neuralyzers
Ayman Akhras
Nathan Luong
Patrick Zhou
Kelly Deng
Reza Jodeiri

January 16, 2025

1 Revision History

Date	Version	Notes
Date 1	1.0	Notes
Date 2	1.1	Notes

2 Symbols, Abbreviations and Acronyms

Symbol	Description
SRS	Software Requirements Specification
AI	Artificial Intelligence
CNN	Convolutional Neural Network
DICOM	Digital Imaging and Communications in Medicine
IVDDs	In Vitro Diagnostic Devices
ML	Machine Learning
PACS	Picture Archiving and Communication System
SaMD	Software as a Medical Device
ROC	Receiver Operating Characteristic Curve
SLA	Service-Level Agreement
FR	Functional Requirement
NFR	Non-Functional Requirement
FSM	Finite State Machine
CXR	Chest X-Ray Project
POC	Proof of Concept
TM	Theoretical Model
AWS	Amazon Web Services
ECS	Elastic Container Service
ECR	Elastic Container Registry

Contents

1	Rev	vision History	i						
2									
3									
4	Not 4.1	tation Data Types from Libraries	1 2						
5		dule Decomposition	3						
6		b Application Server Module	4						
	6.1	Module	4						
	6.2	Uses	4						
	6.3	Syntax	4						
		6.3.1 Exported Constants	4						
	6.4	6.3.2 Exported Access Programs	$\frac{4}{4}$						
	0.4	6.4.1 State Variables	4						
		6.4.2 Environment Variables	4						
		6.4.3 Assumptions	4						
		6.4.4 Access Routine Semantics	5						
		6.4.5 Local Functions	5						
—	TTO	TD Commen Madela	c						
7	7.1	TP Server Module Other Modules the Current Module Uses	6						
	7.1	State Variables	6						
	7.2	Exported Constants and Access Programs	6						
	1.5	7.3.1 Exported Access Programs	6						
		7.3.2 Exported Constants	6						
	7.4	Environment Variables	6						
	7.4 - 7.5	Assumptions	6						
	7.6	Access Routine Semantics	7						
	1.0	7.6.1 startServer()	7						
		7.6.2 stopServer()	7						
		7.6.3 processRequest(request)	7						
	7.7	Local Functions	7						
8	Dig	ease Prediction Server Module	7						
O	8.1	Other Modules the Current Module Uses	7						
	8.2	State Variables	8						
	8.3	Exported Constants and Access Programs	8						

		8.3.1 Exported Access Programs
		8.3.2 Exported Constants
	8.4	Environment Variables
	8.5	Assumptions
	8.6	Access Routine Semantics
		8.6.1 loadModel()
		8.6.2 predictDisease(patientImageData)
	8.7	Local Functions
9	T I a a a	r Authentication Module
9		
	9.1 9.2	Other Modules the Current Module Uses
	9.3	Exported Constants and Access Programs
		9.3.1 Exported Access Programs
	0.4	9.3.2 Exported Constants
	9.4	Environment Variables
	9.5	Assumptions
	9.6	Access Routine Semantics
		9.6.1 sign_up_user(user_email, password, user_detail)
		9.6.2 sign_in_user(user_email, password)
		9.6.3 get_self_user(bearer_token)
		9.6.4 get_user_by_id(uuid)
		9.6.5 update_user_by_id(uuid, updated_user_json)
	9.7	Local Functions
10	Mod	dule NAME HERE!!!
	10.1	Other Modules the Current Module Uses
	10.2	State Variables
		Exported Constants and Access Programs
		10.3.1 Exported Access Programs
		10.3.2 Exported Constants
	10.4	Environment Variables
	10.5	Assumptions
	10.6	Access Routine Semantics
		10.6.1 trackProgression(patientID, data)
		10.6.2
	10.7	Local Functions
11	Doc	tor Profile View Module 14
11		Other Modules the Current Module Uses
		State Variables
		Exported Constants and Access Programs
	11.0	11.3.1 Exported Access Programs

		11.3.2 Exported Constants	14
	11.4	Environment Variables	14
	11.5	Assumptions	15
	11.6	Access Routine Semantics	15
		11.6.1 viewDoctorProfile(doctorID)	15
	11.7	Local Functions	15
12		ent List View Module	15
		Other Modules the Current Module Uses	15
		State Variables	15
	12.3	Exported Constants and Access Programs	16
		12.3.1 Exported Access Programs	16
	40.4	12.3.2 Exported Constants	16
		Environment Variables	16
		Assumptions	16
	12.6	Access Routine Semantics	16
		12.6.1 viewPatientList(doctorID, filters)	16
		12.6.2 applyFilter(filters)	17
		12.6.3 sortList(criteria)	17
	12.7	Local Functions	17
13	Pati	ent Overview Module	17
-0		Other Modules the Current Module Uses	17
		State Variables	17
		Exported Constants and Access Programs	17
	10.0	13.3.1 Exported Access Programs	17
		13.3.2 Exported Constants	18
	13 4	Environment Variables	18
		Assumptions	18
		Access Routine Semantics	18
	10.0	13.6.1 viewPatientOverview(patientID)	18
	13 7	Local Functions	18
	10.1	Econ Turous	10
14	Dise	eases Progression View	18
	14.1	Other Modules the Current Module Uses	18
	14.2	State Variables	19
	14.3	Exported Constants and Access Programs	19
		14.3.1 Exported Access Programs	19
		14.3.2 Exported Constants	19
	14.4	Environment Variables	19
	14.5	Assumptions	19
	14.6	Access Routine Semantics	20
		14.6.1 viewDiseaseProgression(nationtID)	20

1	14.7	Local Functions		 •						. 20
		lical Records List View Module								20
1	15.1	Other Modules the Current Module Uses .								. 20
1	15.2	State Variables								. 20
1	15.3	Exported Constants and Access Programs								. 20
		15.3.1 Exported Access Programs								. 20
		15.3.2 Exported Constants								. 21
1	15.4	Environment Variables								. 21
1	15.5	Assumptions								. 21
1	15.6	Access Routine Semantics								. 21
		15.6.1 viewMedicalRecordsList(patientID)								. 21
1	15.7	Local Functions								
16 I	Med	lical Record View Module								22
1	16.1	Other Modules the Current Module Uses .								. 22
1	16.2	State Variables								. 22
		Exported Constants and Access Programs								
		16.3.1 Exported Access Programs								. 22
		16.3.2 Exported Constants								
1	16.4	Environment Variables								. 22
1	16.5	Assumptions								. 23
		Access Routine Semantics								
		16.6.1 viewMedicalRecord(recordID)								. 23
1	16.7	Local Functions								
17 2	X-R	ay Report View								23
		place holder								. 23
		State Variables								
		Exported Constants and Access Programs								
		17.3.1 Exported Access Programs								
		17.3.2 Exported Constants								. 24
1	17.4	Environment Variables								
		Assumptions								
		Access Routine Semantics								
_		17.6.1 place holder(place holder)								
1	17.7	Local Functions								
		ease Progression Module								2 4
		Other Modules the Current Module Uses .								
		State Variables								
1	10.3	Exported Constants and Access Programs								
		18.3.1 Exported Access Programs								. 25

	18.3.2 Exported Constants						. 2	١
18.4	Environment Variables						. 2	١
18.5	Assumptions							
	Access Routine Semantics							
	18.6.1 loadProgressionModel(modelPath)							
	18.6.2 updateProgressionStage(patientID, stage							
	18.6.3 getProgressionStatus(patientID)							
	18.6.4 resetProgressionModel()							
	18.6.5 forecastProgression(patientID)							
18.7	Local Functions							
	ease Prediction Module						2	
19.1	Other Modules the Current Module Uses							
	State Variables							
19.3	Exported Constants and Access Programs							
	19.3.1 Exported Access Programs						. 2	7
	19.3.2 Exported Constants						. 2	7
19.4	Environment Variables						. 2	7
19.5	Assumptions						. 2	3
19.6	Access Routine Semantics						. 2	3
	19.6.1 loadModel(modelPath)						. 2	3
	19.6.2 predictDisease(patientData)							3
	19.6.3 updateModel(newModel)							3
	19.6.4 getModelStatus()							
	19.6.5 resetModel()							
19.7	Local Functions							96
00 D	B 11 12 12						-	_
	a Persistent Module						2	
	Other Modules the Current Module Uses							
	State Variables							
20.3	Exported Constants and Access Programs							
	20.3.1 Exported Access Programs							
	20.3.2 Exported Constants							
	Environment Variables							
	Assumptions							
20.6	Access Routine Semantics							
	$20.6.1 \text{ get_item_by_id}(\text{uuid}) \dots \dots \dots$							
	20.6.2 create_new_item(item)							
	20.6.3 update_item_by_id(uuid, item)							
	20.6.4 get_s3_upload_presigned_url(file_name)							
	Local Functions							
	Other Modules the Current Module Uses							
20.9	State Variables	 					. 3	2

20.10Exported Constants and Access Programs
20.10.1 Exported Access Programs
20.10.2 Exported Constants
20.11Environment Variables
20.12Assumptions
20.13Access Routine Semantics
$20.13.1 \operatorname{connectDB}(\operatorname{dbPath}) \dots 33$
$20.13.2 \text{saveData}(\text{data}) \dots 33.2 \text{saveData}(\text{data}) \dots 33.2 \text{saveData}(\text{data}) \dots \dots$
$20.13.3 \text{loadData(query)} \dots \dots \dots 33.3 \text{loadData(query)} \dots \dots \dots \dots \dots \dots \dots \dots \dots $
$20.13.4\mathrm{updateData(dataID,newData)}$
$20.13.5 \text{deleteData}(\text{dataID}) \dots 34$
$20.13.6\mathrm{backupDatabase}(\mathrm{backupPath})$
$20.13.7 \mathrm{restoreDatabase}(\mathrm{backupPath}) \dots 34$
20.14Local Functions

3 Introduction

The following document details the Module Interface Specifications for the application. Complementary documents include the Module Guide.

The full documentation and implementation can be found at CXR-Capstone.

4 Notation

The structure of the MIS for modules comes from Hoffman and Strooper (1995). The mathematical notation follows standard conventions for sets, sequences, and functions. The notation includes symbols for ML model operations and medical data processing.

Data Type	Notation	Description
character	char	a single symbol or digit
integer	\mathbb{Z}	a number without a fractional component in $(-\infty, \infty)$
natural number	N	a number without a fractional component in $[1, \infty)$
real	\mathbb{R}	any number in $(-\infty, \infty)$
tensor	\mathbb{T}^n	an n-dimensional array of numerical values used for ML computations
probability	\mathbb{P}	a real number in $[0, 1]$ representing likelihood
matrix	$\mathbb{M}_{m,n}$	a 2D array of size $m \times n$ containing numerical values
binary	\mathbb{B}	boolean values true, false

4.1 Data Types from Libraries

Data Type	Notation	Description
DICOM	DICOM	Digital Imaging and Communications in Medicine format, standard for medical imaging storage and transmission
PyTorch Tensor	torch.Tensor	Multi-dimensional matrix containing elements of a single data type, optimized for GPU operations
NumPy Array	np.ndarray	Multi-dimensional array for scientific computing and image processing
JWT	JWT	JSON Web Token for secure authentication and information transmission
HTTP Request	HTTPRequest	Object containing HTTP request information including headers, body, and method
HTTP Response	HTTPResponse	Object containing HTTP response information including status code, headers, and body
DataFrame	pd.DataFrame	2-dimensional labeled data structure for patient and medical records
Base64	Base64	Binary-to-text encoding scheme for transmitting binary image data
YAML	YAML	Human-readable data serialization format for configuration files

Table 1: Data types from libraries

5 Module Decomposition

The following table is taken directly from the Module Guide document for this project.

Level 1	Level 2
Hardware-Hiding Module	Web Application Server
	HTTP Server
	Disease Prediction Server
	Disease Progression Server
	User Authentication
	Patients List View
	Patient Overview View
Behaviour-Hiding Module	Disease Progression View
	Medical Records List View
	X-Ray Report View
	Disease Progression Model
Software Decision Module	Disease Prediction Model
	Data Persistent

Table 2: Module Hierarchy

6 Web Application Server Module

6.1 Module

Web Application Server

6.2 Uses

N/A

6.3 Syntax

6.3.1 Exported Constants

N/A

6.3.2 Exported Access Programs

Name	In	Out	Exceptions
handleRequest	HTTP request	HTTP response	Invalid Request Exception

6.4 Semantics

6.4.1 State Variables

- sessionData: Stores current session information.
- activeUsers: Keeps track of currently active users.

6.4.2 Environment Variables

- serverPort: The port on which the server listens for incoming connections.
- hostAddress: The server's host address.

6.4.3 Assumptions

- Assumes the HTTP Server Module (M2) is properly configured and running.
- Assumes valid HTTP requests are received.

6.4.4 Access Routine Semantics

handleRequest(request)

- transition: Processes the incoming HTTP request and routes it to the appropriate view module.
- output: Returns the HTTP response based on the request.

6.4.5 Local Functions

- parseRequest(request): Parses the incoming HTTP request to extract necessary information.
- generateResponse(data): Constructs an HTTP response based on the processed data.
- authenticateUser(credentials): Verifies the user's credentials before processing the request.

7 HTTP Server Module

7.1 Other Modules the Current Module Uses

• N/A

7.2 State Variables

• requestQueue: Holds incoming HTTP requests until they are processed.

• responseQueue: Holds outgoing HTTP responses that need to be sent back to clients.

7.3 Exported Constants and Access Programs

7.3.1 Exported Access Programs

Name	In	Out	Exceptions
startServer	None	None	ServerStartException
stopServer	None	None	ServerStopException
processRequest	HTTP request	HTTP response	InvalidRequestException

7.3.2 Exported Constants

• **SERVER_PORT**: 8080

• MAX_CONNECTIONS: 100

7.4 Environment Variables

• serverPort: The port on which the server listens for incoming HTTP connections.

• maxConnections: Maximum number of simultaneous connections the server can handle.

7.5 Assumptions

- Assumes the Web Application Server Module (M1) is properly configured and running.
- Assumes incoming HTTP requests are formatted correctly.

7.6 Access Routine Semantics

7.6.1 startServer()

- Transition: Starts the HTTP server, initializes necessary resources, and begins listening for incoming requests.
- Output: No output, but may throw a ServerStartException if the server cannot be started.

7.6.2 stopServer()

- Transition: Stops the HTTP server, gracefully shuts down connections.
- Output: No output, but may throw a ServerStopException if the server cannot be stopped.

7.6.3 processRequest(request)

- **Transition**: Takes an incoming HTTP request and processes it, routing it to the appropriate server module or view module.
- Output: Returns an HTTP response based on the processed request.

7.7 Local Functions

- parseRequest(): Parses the incoming HTTP request to extract necessary information such as headers and parameters.
- generateResponse(): Constructs an HTTP response based on the processed data from the request.
- handleError(): Handles errors that arise during request processing and generates appropriate error responses.

8 Disease Prediction Server Module

8.1 Other Modules the Current Module Uses

- M1: Web Application Server Module
- M2: HTTP Server Module
- M14: Disease Prediction Model

8.2 State Variables

- model: The pre-trained model from torchxrayvision used for predicting lung diseases from X-ray images.
- modelAccuracy: Tracks the accuracy of the current model after training and validation.
- **predictionThreshold**: A constant threshold to determine the classification outcome (e.g., disease presence).
- patientImageData: Holds the chest X-ray image data used for prediction.

8.3 Exported Constants and Access Programs

8.3.1 Exported Access Programs

Name	In	Out	Exceptions
loadModel	None	Loaded model	ModelLoadException
predictDisease	X-ray image data	Disease prediction	InvalidImageException

8.3.2 Exported Constants

- PREDICTION_THRESHOLD: 0.75 (threshold for classification of disease presence)
- MODEL_PATH: Path to the pre-trained model (e.g., ./models/chest_xray_model.pth)
- MAX_PREDICTIONS: 1000 (maximum number of predictions to handle concurrently)

8.4 Environment Variables

- modelPath: The path where the torchxrayvision pre-trained model is saved or loaded from.
- **predictionEndpoint**: The endpoint for making predictions using chest X-ray images.

8.5 Assumptions

- Assumes the pre-trained torchxrayvision model is available and compatible with the data provided.
- Assumes valid X-ray image data is available for predictions.
- Assumes the Web Application Server Module (M1) and HTTP Server Module (M2) are properly configured and running.

8.6 Access Routine Semantics

8.6.1 loadModel()

• **Transition**: Loads the pre-trained disease prediction model from the specified path after image is uploaded using torchxrayvision.

8.6.2 predictDisease(patientImageData)

- **Transition**: Uses the loaded model to make predictions based on the provided X-ray image data.
- Output: Returns the disease prediction (e.g., probability of a disease being present) or throws an InvalidImageException if the image is invalid.

8.7 Local Functions

- loadModel(): Loads the pre-trained model from disk or cloud storage using torchxrayvision's functionality.
- evaluateModel(): Evaluates the model's performance with a test dataset to calculate metrics like accuracy and sensitivity.
- **preprocessImage()**: Preprocesses incoming X-ray image data to fit the model's input requirements (e.g., resizing, normalization).
- postprocessPrediction(): Processes the raw output from the model (e.g., probabilities) into a human-readable format (e.g., disease labels).

9 User Authentication Module

9.1 Other Modules the Current Module Uses

• M15: Data Persistence Module

9.2 State Variables

- cognito_idp_client: AWS Cognito client, initialized via the AWS boto3 library.
- user_pool_id: Identifier of the database at which AWS Cognito stores user's information such as username, password, first name, etc.
- client_id: Identifier of the User Authentication Module to authenticate with AWS Cognito (via OAUTH's client credential grant).
- **client_secret**: Password of the User Authentication Module, used in conjunction with client_id to authenticate with AWS Cognito.

9.3 Exported Constants and Access Programs

9.3.1 Exported Access Programs

Name	In	Out	Exceptions
sign_up_user	user_email, password, user_detail	user	Unauthorized, UserAlreadyEx
sign_in_user	user_email, password	bearer_token	Unauthorized, IncorrectCred
get_self_user	bearer_token	user	Unauthorized, InternalServe
get_user_by_id	uuid	user	Unauthorized, UserNotFoundE
update_user_by_id	uuid, updated_user_json	user	Unauthorized, UserNotFoundE

9.3.2 Exported Constants

This module does not export any constants

9.4 Environment Variables

- AWS_ACCESS_KEY_ID: Access key to authenticate backend with AWS
- AWS_SECRET_ACCESS_KEY: Secret key to authenticate backend with AWS
- AWS_REGION: AWS Region of the backend infrastructure

9.5 Assumptions

- AWS Cognito returns consistent json results with user fields, since neuralanalyzer backend expect certain fields from user object (custom:firstName, custom:organization, etc.)
- AWS Cognito always stays online, since without this service, users wont be able to login into the application.

9.6 Access Routine Semantics

9.6.1 sign_up_user(user_email, password, user_detail)

- Inputs:
 - user_email: Email address of the user.
 - password: Password for the user's account.
 - user_detail: Additional user details, such as first name and organization, formatted as Python Dictionary.
- Outputs: A user object containing the created user's details.

• Preconditions:

- The provided email is valid and not already registered.
- Password meets AWS Cognito's security requirements.

• Postconditions:

- The user is successfully created and stored in the Cognito user pool.

9.6.2 sign_in_user(user_email, password)

• Inputs:

- **user_email**: Email address of the user.
- **password**: Password for the user's account.
- Outputs: A bearer token for authenticated access to other services.

• Preconditions:

- The user is registered in the Cognito user pool.
- The provided credentials are correct.

• Postconditions:

- The user is successfully authenticated.
- A bearer token is generated and returned for future authenticated requests.

9.6.3 get_self_user(bearer_token)

• Inputs:

- bearer_token: User's Bearer token.
- Outputs: A user object containing the created user's details.

• Preconditions:

- The bearer_token must be valid and not expired.

• Postconditions:

- The user is successfully retrieved from AWS Cognito.

9.6.4 get_user_by_id(uuid)

• Inputs:

- **uuid**: The unique identifier of the user (string)
- Outputs: A user object containing the user's details.

• Preconditions:

- User's UUID must be valid, and attached to a registered user.

• Postconditions:

- The user is successfully retrieved from AWS Cognito.

9.6.5 update_user_by_id(uuid, updated_user_json)

• Inputs:

- **uuid**: The unique identifier of the user (string)
- updated_user_json: Python Dictionary object that contains the field which the caller wants to update the user with.
- Outputs: A user object containing the updated user's details.

• Preconditions:

- User's UUID must be valid, and attached to a registered user.
- User's updated fields must belong to list of fields users are allowed to edit (ie: first and last name, role, organization, etc.)

• Postconditions:

 User with UUID must have the fields in updated_user_json updated on AWS Cognito.

9.7 Local Functions

- _parse_user_json_from_cognito_user: Convert user object returned from AWS Cognito, to a more application-friendly internal user object.
- _parse_cognito_attr_user_from_user_json: Convert the internal user object to the user data format that AWS Cognito requires.
- _get_username_from_email: Converting user's email into their username. (ie: nathan@gmail.com will be converted to nathan_gmail)
- _get_user_by_email: Fetching user from AWS Cognito via their email.

10 Module NAME HERE!!!

10.1 Other Modules the Current Module Uses

• fill this

10.2 State Variables

• Title: fill this

10.3 Exported Constants and Access Programs

10.3.1 Exported Access Programs

Name	In	Out	Exceptions
fill 1	fill 2	fill 3	fill 4
fill 1	fill 2	fill 3	fill 4

10.3.2 Exported Constants

• Title: file this

10.4 Environment Variables

• fill this

10.5 Assumptions

• fill this

10.6 Access Routine Semantics

10.6.1 trackProgression(patientID, data)

• fill this

10.6.2

• fill this

10.7 Local Functions

• fill this

11 Doctor Profile View Module

11.1 Other Modules the Current Module Uses

• M1: Web Application Server Module

• M2: HTTP Server Module

• M5: User Authentication Module

• M7: Patients Overview View

• M10: Medical Records List View

• M15: Data Persistence Module

11.2 State Variables

• **doctorProfile**: Stores the user information of the currently logged-in doctor, including name, specialty, contact details, and credentials.

11.3 Exported Constants and Access Programs

11.3.1 Exported Access Programs

Name	In	Out	Exceptions
getDoctorProfile	userID	user Profile	ProfileNotFoundException
getDoctorProfile	userID	user Profile	NotDoctorException

11.3.2 Exported Constants

- PROFILE_UPDATE_INTERVAL: 24 hours (time interval for automatic profile updates)
- **DEFAULT_PROFILE_PICTURE**: "default_doctor.png" (default profile picture for doctors without a custom one)
- **DEFAULT_SPECIALTY**: "Radiologist" (default specialty for doctors without a specified specialty)

11.4 Environment Variables

- profileDataPath: Path to the data source containing doctor profile information.
- authenticationService: URL or endpoint of the service used for user authentication.

11.5 Assumptions

- Assumes that the User Authentication Module (M5) ensures only authorized doctors can view their profiles.
- Assumes the profile data is accurately stored and updated in the Data Persistence Module (M15).

11.6 Access Routine Semantics

11.6.1 viewDoctorProfile(doctorID)

- **Transition**: Retrieves the profile information of the specified doctor.
- Output: Returns the doctor profile details or throws a ProfileNotFoundException if the profile cannot be found.

11.7 Local Functions

- fetchProfileData(doctorID): Retrieves profile data from the data source.
- updateProfilePicture(doctorID, picturePath): Updates the profile picture of the doctor.
- logProfileAccess(doctorID): Logs each time a doctor accesses their profile for auditing purposes.

12 Patient List View Module

12.1 Other Modules the Current Module Uses

- M1: Web Application Server Module
- M2: HTTP Server Module
- M5: User Authentication
- M7: Patient Overview View
- M15: Data Persistent

12.2 State Variables

- patientList: Stores a list of patients assigned to a doctor, including their basic information such as name, age, and medical condition.
- searchFilters: Stores the current filters applied to the patient list for sorting and searching purposes.

12.3 Exported Constants and Access Programs

12.3.1 Exported Access Programs

Name	In	Out	Exceptions
viewPatientList	Doctor ID, filters	List of patient details	PatientListNotFoundException
applyFilter	Filter parameters	Filtered patient list	InvalidFilterException
sortList	Sorting criteria	Sorted patient list	None

12.3.2 Exported Constants

- MAX_PATIENTS_PER_PAGE: 20 (maximum number of patients displayed per page in the list)
- **DEFAULT_SORT_ORDER**: "alphabetical" (default order in which patients are listed)

12.4 Environment Variables

- patientDataPath: Path to the data source containing patient records.
- authenticationService: URL or endpoint of the service used for user authentication.

12.5 Assumptions

- Assumes that the User Authentication Module (M5) is responsible for verifying that the user is a doctor with access to the list.
- Assumes the patient data is up-to-date and synchronized with the Data Persistence Module (M15).

12.6 Access Routine Semantics

12.6.1 viewPatientList(doctorID, filters)

- **Transition**: Retrieves a list of patients associated with the given doctor, applying the specified filters.
- Output: Returns the list of patient details or throws a PatientListNotFoundException if no patients are found.

12.6.2 applyFilter(filters)

- Transition: Applies the given filters to the current patient list.
- Output: Returns the filtered patient list or throws an InvalidFilterException if the filters are invalid.

12.6.3 sortList(criteria)

- Transition: Sorts the current patient list based on the provided criteria.
- Output: Returns the sorted patient list.

12.7 Local Functions

- filterPatients(filters): Filters the patient list according to the specified parameters.
- sortPatients(criteria): Sorts the patient list based on the provided criteria.
- logListAccess(doctorID): Logs each time a doctor accesses the patient list for auditing purposes.

13 Patient Overview Module

13.1 Other Modules the Current Module Uses

- M1: Web Application Server Module
- M2: HTTP Server Module
- M11: User Authentication Module
- M15: Data Persistence Module

13.2 State Variables

• patientOverview: Stores the summary information of a selected patient, including personal details, recent medical history, and current treatment plans.

13.3 Exported Constants and Access Programs

13.3.1 Exported Access Programs

Name	In	Out	Exceptions
getPatient	Patient ID	Patient overview details	OverviewNotFoundException

13.3.2 Exported Constants

• OVERVIEW_REFRESH_INTERVAL: 15 minutes (interval for refreshing the patient overview data)

13.4 Environment Variables

- overview DataPath: Path to the data source containing patient overview information.
- authenticationService: URL or endpoint of the service used for user authentication.

13.5 Assumptions

- Assumes the User Authentication Module (M5) ensures that only authorized users can view the patient overview.
- Assumes the overview data is accurate and synchronized with the Data Persistence Module (M15).

13.6 Access Routine Semantics

13.6.1 viewPatientOverview(patientID)

- Transition: Retrieves the overview information of the specified patient.
- Output: Returns the patient overview details or throws an OverviewNotFoundException if the overview cannot be found.

13.7 Local Functions

- fetchOverviewData(patientID): Retrieves overview data from the data source.
- **generateOverviewSummary(patientID)**: Creates a summary of the patient's current status.
- logOverviewAccess(patientID): Logs each time a patient overview is accessed for auditing purposes.

14 Diseases Progression View

14.1 Other Modules the Current Module Uses

- M1: Web Application Server Module
- M2: HTTP Server Module

• M11: User Authentication Module

• M12: X-Ray Report View

• M15: Data Persistence Module

14.2 State Variables

• Title: fill this

14.3 Exported Constants and Access Programs

14.3.1 Exported Access Programs

Name	In	Out	Exceptions
viewProgression	Patient ID	Progression data	ProgressionNotFoundException

14.3.2 Exported Constants

- PROGRESSION_UPDATE_INTERVAL: 12 hours (time interval for updating disease progression data)
- **DEFAULT_CHART_TEMPLATE**: "default_progression_chart.html" (default template for displaying disease progression charts)

14.4 Environment Variables

- **progressionDataPath**: Path to the data source containing disease progression information.
- **chartRenderingService**: URL or endpoint of the service used for rendering progression charts.

14.5 Assumptions

- Assumes the User Authentication Module (M11) ensures that only authorized users can view the disease progression.
- Assumes the data is regularly updated and maintained in the Data Persistence Module (M15).

14.6 Access Routine Semantics

14.6.1 viewDiseaseProgression(patientID)

- **Transition**: Retrieves the disease progression details for the specified patient.
- Output: Returns the progression details or throws a ProgressionNotFoundException if the data cannot be found.

14.7 Local Functions

- **fetchProgressionData(patientID)**: Retrieves disease progression data from the data source.
- **generateProgressionChart(patientID)**: Generates a visual chart of the disease progression.
- logProgressionAccess(patientID): Logs access to a patient's disease progression data for auditing purposes.

15 Medical Records List View Module

15.1 Other Modules the Current Module Uses

- M1: Web Application Server Module
- M2: HTTP Server Module
- M5: User Authentication Module
- M11: Medical Record View Module
- M15: Data Persistence Module

15.2 State Variables

• medicalRecordsList: Contains a list of all medical records associated with a specific patient, including dates, record types, and summaries.

15.3 Exported Constants and Access Programs

15.3.1 Exported Access Programs

Name	In	Out	Exceptions
viewMedicalRecordsList	Patient ID	List of medical records	RecordsNotFoundException

15.3.2 Exported Constants

- RECORDS_REFRESH_INTERVAL: 10 minutes (interval for refreshing the list of medical records)
- **DEFAULT_LIST_TEMPLATE**: "default_records_list_template.html" (default template for displaying the list of medical records)

15.4 Environment Variables

- recordsDataPath: Path to the data source containing the patient's medical records.
- recordSummaryService: URL or endpoint of the service used for summarizing medical records.

15.5 Assumptions

- Assumes the User Authentication Module (M5) ensures that only authorized users can view the medical records.
- Assumes the records data is accurate and up-to-date in the Data Persistence Module (M15).

15.6 Access Routine Semantics

15.6.1 viewMedicalRecordsList(patientID)

- Transition: Retrieves the list of medical records for the specified patient.
- Output: Returns the list of medical records or throws a RecordsNotFoundException if no records are found.

15.7 Local Functions

- **fetchMedicalRecords(patientID)**: Retrieves the list of medical records from the data source.
- generateRecordsSummary(patientID): Creates summaries for each record in the list.
- logRecordsAccess(patientID): Logs access to a patient's medical records list for auditing purposes.

16 Medical Record View Module

16.1 Other Modules the Current Module Uses

• M1: Web Application Server Module

• M2: HTTP Server Module

• M11: User Authentication Module

• M12: X-Ray Report View

• M15: Data Persistence Module

16.2 State Variables

• medicalRecordDetails: Stores detailed information about a specific medical record, including the date, diagnosis, treatment, and doctor notes.

16.3 Exported Constants and Access Programs

16.3.1 Exported Access Programs

Name	In	Out	Exceptions
viewMedicalRecord	Record ID	Medical record details	RecordNotFoundException

16.3.2 Exported Constants

- **DEFAULT_RECORD_TEMPLATE**: "default_record_template.html" (default template for displaying a medical record)
- MAX_RECORD_DISPLAY_LENGTH: 5000 characters (maximum length for displaying record content)

16.4 Environment Variables

- recordDataPath: Path to the data source containing the specific medical record.
- recordDetailService: URL or endpoint of the service used for fetching detailed medical record data.

16.5 Assumptions

- Assumes the User Authentication Module (M5) ensures that only authorized users can view the detailed medical record.
- Assumes the data for medical records is accurately maintained and readily accessible through the Data Persistence Module (M15).

16.6 Access Routine Semantics

16.6.1 viewMedicalRecord(recordID)

- Transition: Retrieves the details of the specific medical record identified by recordID.
- Output: Returns the medical record details or throws a RecordNotFoundException if the record cannot be found.

16.7 Local Functions

- **fetchMedicalRecord(recordID)**: Retrieves the details of the specified medical record from the data source.
- formatRecordDetails(recordID): Formats the medical record details for display.
- logRecordAccess(recordID): Logs access to the specific medical record for auditing purposes.

17 X-Ray Report View

17.1 place holder

- M1: Web Application Server Module
- M2: HTTP Server Module
- M15: User Authentication Module
- M13: Disease Prediction Server Module
- M15: Data Persistence Server Module

17.2 State Variables

• place holder: place holder

17.3 Exported Constants and Access Programs

17.3.1 Exported Access Programs

Name	In	Out	Exceptions
place holder	place holder	place holder	place holder

17.3.2 Exported Constants

• place holder: "place holder"

17.4 Environment Variables

• place holder: place holder

17.5 Assumptions

- place holder
- place holder

17.6 Access Routine Semantics

17.6.1 place holder(place holder)

• place holder: place holder

17.7 Local Functions

• place holder(place holder): place holder

18 Disease Progression Module

18.1 Other Modules the Current Module Uses

- M1: Web Application Server Module
- M2: HTTP Server Module
- M5: User Authentication Module
- M15: Data Persistence Module

18.2 State Variables

currentProgressionModel: Model – The currently loaded disease progression model.

progressionStatus: ProgressionStatus - Status of the progression model (e.g., loaded, un-

loaded, error).

progressionCache: Cache – Caches recent progression results for quick access.

18.3 Exported Constants and Access Programs

18.3.1 Exported Access Programs

Name	In	Out	Exception
loadProgressionModel	modelPath : String	-	ModelLoad
updateProgressionStage	patientID : String, stageData : StageData	-	InvalidSt
getProgressionStatus	patientID : String	ProgressionStatus	PatientNo
resetProgressionModel	-	-	-
forecastProgression	patientID : String	ForecastResult	ForecastE

18.3.2 Exported Constants

- **DEFAULT_PROGRESSION_MODEL**: "progression_model.pkl" Path to the default progression model.
- MAX_STAGES: 10 Maximum number of disease progression stages.

18.4 Environment Variables

- PROGRESSION_MODEL_PATH: Path to the directory containing progression models.
- CACHE_SIZE: Maximum number of progression forecasts to cache.

18.5 Assumptions

- Assumes that the Data Persistence Module (M15) provides reliable access to patient data.
- Assumes that the User Authentication Module (M5) ensures that only authorized users can update and view disease progression.
- Assumes that progression models are regularly updated and maintained for accuracy.

18.6 Access Routine Semantics

18.6.1 loadProgressionModel(modelPath)

- Transition: Loads the disease progression model from the specified modelPath.
- Exception: Throws ModelLoadException if the model cannot be loaded.

18.6.2 updateProgressionStage(patientID, stageData)

- **Transition**: Updates the disease progression stage for the patient identified by **patientID** with the provided **stageData**.
- Exception: Throws InvalidStageException if stageData is invalid.

18.6.3 getProgressionStatus(patientID)

- **Transition**: Retrieves the current disease progression status for the patient identified by patientID.
- Output: Returns a ProgressionStatus object.
- Exception: Throws PatientNotFoundException if the patient does not exist.

18.6.4 resetProgressionModel()

- Transition: Resets the disease progression model to the default state.
- Exception: None.

18.6.5 forecastProgression(patientID)

- **Transition**: Generates a forecast of the disease progression for the patient identified by patientID.
- Output: Returns a ForecastResult containing predicted progression data.
- Exception: Throws ForecastException if the forecast cannot be generated.

18.7 Local Functions

- validateStageData(stageData): Validates the structure and completeness of stageData.
- loadDefaultProgressionModel(): Loads the default progression model.
- cacheProgressionResult(result): Stores the ForecastResult in the progressionCache.
- clearProgressionCache(): Clears all entries from the progressionCache.

19 Disease Prediction Module

19.1 Other Modules the Current Module Uses

• M1: Web Application Server Module

• M2: HTTP Server Module

• M5: User Authentication Module

• M15: Data Persistence Module

19.2 State Variables

• currentModel: Model and The currently loaded disease prediction model.

• modelStatus: ModelStatus and Status of the prediction model (e.g., loaded, unloaded, error).

• predictionCache: Cache and Caches recent prediction results for quick access.

19.3 Exported Constants and Access Programs

19.3.1 Exported Access Programs

Name	In	Out	Exceptions
loadModel	modelPath : String	-	ModelLoadException
predictDisease	patientData : PatientData	PredictionResult	InvalidInputException
updateModel	newModel : Model	-	ModelUpdateException
getModelStatus	-	ModelStatus	-
resetModel	-	-	-

19.3.2 Exported Constants

- **DEFAULT_PREDICTION_MODEL**: "default_model.pkl" (path to the default prediction model)
- MAX_INPUT_SIZE: 1024 (maximum size for input data)

19.4 Environment Variables

- MODEL_PATH: Path to the directory containing prediction models.
- CACHE_SIZE: Maximum number of predictions to cache.

19.5 Assumptions

- Assumes that the Data Persistence Module (M15) provides reliable access to patient data.
- Assumes that the User Authentication Module (M11) ensures that only authorized users can perform predictions.
- Assumes that prediction models are compatible with the system's architecture and dependencies.

19.6 Access Routine Semantics

19.6.1 loadModel(modelPath)

- Transition: Loads the disease prediction model from the specified modelPath.
- Exception: Throws ModelLoadException if the model cannot be loaded.

19.6.2 predictDisease(patientData)

- **Transition**: Processes **patientData** using the current prediction model to generate a prediction.
- Output: Returns a PredictionResult containing the predicted disease information.
- Exception: Throws InvalidInputException if patientData is malformed or incomplete.

19.6.3 updateModel(newModel)

- Transition: Updates the current prediction model with newModel.
- Exception: Throws ModelUpdateException if the update fails.

19.6.4 getModelStatus()

- Output: Returns the current status of the prediction model (ModelStatus).
- Exception: None.

19.6.5 resetModel()

- Transition: Resets the prediction model to the default state.
- Exception: None.

19.7 Local Functions

- validatePatientData(patientData): Validates the structure and completeness of patientData.
- loadDefaultModel(): Loads the default prediction model.
- cachePrediction(result): Stores the PredictionResult in the predictionCache.
- clearCache(): Clears all entries from the predictionCache.

20 Data Persistent Module

20.1 Other Modules the Current Module Uses

• None

20.2 State Variables

- **dynamodb**: AWS DynamoDB client (managed NoSQL service), initialized via aws boto3 library.
- table: Name of the DynamoDB table.
- **s3_client**: AWS S3 client (managed blob storage service), initialized via aws boto3 library.

20.3 Exported Constants and Access Programs

20.3.1 Exported Access Programs

Name	In	Out	Exceptions
get_item_by_id	uuid	item	InternalServerError, ClientError
create_new_item	item	uuid, item	InternalServerError, ClientError
update_item_by_id	uuid, item	item	InternalServerError, ClientError
get_s3_upload_presigned_url	file_name	presigned_url	InternalServerError, ClientError

20.3.2 Exported Constants

This module does not export any constants

20.4 Environment Variables

- AWS_ACCESS_KEY_ID: Access key to authenticate backend with AWS
- AWS_SECRET_ACCESS_KEY: Secret key to authenticate backend with AWS
- AWS_REGION: AWS Region of the backend infrastructure

20.5 Assumptions

- AWS DynamoDB and AWS S3 returns consistent json results and error codes.
- AWS DynamoDB and AWS S3 always stays online, without these services, users wont be able to create, view, update prescriptions and records.

20.6 Access Routine Semantics

20.6.1 get_item_by_id(uuid)

- Inputs:
 - **uuid**: Unique identifier of the item to retrieve.
- Outputs:
 - The item object corresponding to the given 'uuid'.
- Preconditions:
 - The 'uuid' exists in the DynamoDB table.
- Postconditions:
 - The item object is successfully retrieved from the database.

20.6.2 create_new_item(item)

- Inputs:
 - item: A Python Dictionary object containing the details of the item to be created.
- Outputs:
 - **uuid**: The unique identifier of the created item.
 - item: The created item object.
- Preconditions:
 - The 'item' dictionary is well-formed.

• Postconditions:

- The new item is successfully added to the DynamoDB table.

20.6.3 update_item_by_id(uuid, item)

• Inputs:

- **uuid**: The unique identifier of the item to be updated.
- item: A Python Dictionary containing the updated fields for the item.

• Outputs:

- The updated item object.

• Preconditions:

- The 'uuid' corresponds to an existing item in the DynamoDB table.
- The 'item' Dictionary contains valid fields for the update.

• Postconditions:

- The item is successfully updated in the database, and returned to user.

20.6.4 get_s3_upload_presigned_url(file_name)

• Inputs:

- file_name: The name of the file to be uploaded to S3 (string).

• Outputs:

 presigned_url: A pre-signed URL allowing authorized upload of the specified file to AWS S3.

• Preconditions:

- The 'file_name' is a valid, non-empty string
- The S3 bucket configured in the s3_client is accessible and writable by the authenticated user.

• Postconditions:

- A valid pre-signed URL for uploading the specified file is generated and returned.
- The generated pre-signed URL is only valid for 1 minute.

20.7 Local Functions

• None

20.8 Other Modules the Current Module Uses

• M1: Web Application Server Module

• M2: HTTP Server Module

20.9 State Variables

dbConnection: DBConnection – Represents the active connection to the database.

connectionStatus: ConnectionStatus – Indicates the status of the database connection.

retryCount: int – Tracks the number of retry attempts for failed operations.

20.10 Exported Constants and Access Programs

20.10.1 Exported Access Programs

Name	In	Out	Exceptions
connectDB	dbPath : String	ConnectionStatus	ConnectionExceptio
saveData	data : DataObject	-	SaveException
loadData	query : Query	DataObject	LoadException
updateData	dataID : String, newData : DataObject	-	UpdateException
deleteData	dataID : String	-	DeleteException
backupDatabase	backupPath : String	-	BackupException
restoreDatabase	backupPath : String	_	RestoreException

20.10.2 Exported Constants

- DEFAULT_DB_PATH: "/var/data/persistence.db" Default database file path.
- MAX_RETRIES: 5 Maximum number of retry attempts for database operations.

20.11 Environment Variables

- **DB_PATH**: Path to the primary database file.
- BACKUP_PATH: Path to store database backups.
- RETRY_LIMIT: Maximum number of retry attempts for database operations.

20.12 Assumptions

- Assumes that the database server is accessible and properly configured.
- Assumes that sufficient storage is available for data persistence and backups.
- Assumes that the User Authentication Module (M11) handles access control for database operations.

20.13 Access Routine Semantics

20.13.1 connectDB(dbPath)

- Transition: Establishes a connection to the database located at dbPath.
- Output: Returns the current ConnectionStatus.
- Exception: Throws ConnectionException if the connection fails.

20.13.2 saveData(data)

- Transition: Saves the provided data object to the database.
- Exception: Throws SaveException if the data cannot be saved.

20.13.3 loadData(query)

- Transition: Executes the query to retrieve data from the database.
- Output: Returns the resulting DataObject.
- Exception: Throws LoadException if the data cannot be retrieved.

20.13.4 updateData(dataID, newData)

- Transition: Updates the data entry identified by dataID with newData.
- Exception: Throws UpdateException if the update fails.

20.13.5 deleteData(dataID)

- Transition: Deletes the data entry identified by dataID from the database.
- Exception: Throws DeleteException if the deletion fails.

20.13.6 backupDatabase(backupPath)

- Transition: Creates a backup of the database at the specified backupPath.
- Exception: Throws BackupException if the backup process fails.

20.13.7 restoreDatabase(backupPath)

- Transition: Restores the database from the backup located at backupPath.
- Exception: Throws RestoreException if the restoration process fails.

20.14 Local Functions

- validateData(data): Validates the structure and integrity of data.
- retryOperation(operation): Attempts to retry a failed operation up to MAX_RETRIES.
- logDatabaseActivity(activity): Logs database operations for auditing and debugging purposes.
- initializeConnection(): Initializes the database connection using environment variables.

References

Daniel M. Hoffman and Paul A. Strooper. Software Design, Automated Testing, and Maintenance: A Practical Approach. International Thomson Computer Press, New York, NY, USA, 1995. URL http://citeseer.ist.psu.edu/428727.html.

Appendix — Reflection

The information in this section will be used to evaluate the team members on the graduate attribute of Problem Analysis and Design.

- 1. What went well while writing this deliverable?
- 2. What pain points did you experience during this deliverable, and how did you resolve them?
- 3. Which of your design decisions stemmed from speaking to your client(s) or a proxy (e.g. your peers, stakeholders, potential users)? For those that were not, why, and where did they come from?
- 4. While creating the design doc, what parts of your other documents (e.g. requirements, hazard analysis, etc), it any, needed to be changed, and why?
- 5. What are the limitations of your solution? Put another way, given unlimited resources, what could you do to make the project better? (LO_ProbSolutions)
- 6. Give a brief overview of other design solutions you considered. What are the benefits and tradeoffs of those other designs compared with the chosen design? From all the potential options, why did you select the documented design? (LO_Explores)

- 1. When writing this deliverable, several things went really well that made the whole process smoother and more organized. Our existing documentation, like the detailed list of functional and non-functional requirements and the Module Interface Specification (MIS) we created, provided a solid roadmap for how the system should behave. This clarity made it easier to develop the design blueprint and the specifications for each module because we already knew what each part needed to do and how they should interact. Another big plus was the team's early focus on a modular architecture, such as separating the machine learning (ML) services from the core backend. This approach kept our design documents organized and manageable, allowing us to develop, test, and maintain each module independently without getting overwhelmed. Additionally, having the core implementation set up during the Proof of Concept (POC) phase gave us a great direction for developing the specifications for our modules. It provided a tested foundation to build on, reducing uncertainties and making our specification development more precise and targeted.
- 2. One of the pain points when writing this document is the ambiguities in some of our existing documents which made it challenging for us to finalize module boundaries. To resolve this, we had couple discussions as a team to refine unclear requirements and ensure alignment. We also experienced a hard time determining the design architecture of our project between MVC and layered architecture. MVC is primarily used for organizing applications with user interfaces, especially in web and desktop applications, whereas layered architecture is highly scalable and maintainable and is easy to replace or upgrade individual layers. At the end we decided to go with layered architecture since when building larger-scale, complex systems that need clear abstractions and the system will integrate with multiple services or technologies. When multiple team members were contributing to the same document, we had a hard time managing version control. To resolve this, we implemented a review-and-approval workflow, which each commit needs at least one another team member's approval before merging in.
- 3. After speaking with our client, we decided to prioritize user-friendly interface which was based on client (professor) feedback as it needs to be intuitive and easy to navigate, which also emphasizes maintainability of the project. We decided to prioritize interpretability and transparency of AI-generated report (explainable AI) which comes from the discussions with stakeholder emphasizing trust in AI results. This would help doctors/radiologists to better understand the prediction results. As for non-client-informed decisions, we decided to implement a disease progression feature that is separated from the prediction model. It traces the patient disease progression over a certain period, allowing comparing different X-rays and reports to give a more straightforward understanding of patient condition overtime. By separating two models, they work independently and will still function if the other one goes down, allowing higher fault tolerance.
- 4. While creating this design doc, as we investigated the functionality and operational

flow of the system, we recognized gaps in our initial architecture that required additional components to ensure a logical and coherent flow. We found out that the Software Requirements Specification and hazard analysis were incomplete and needed to be revisited to accommodate these changes. For SRS we might need to adjust the specification for inputs taken into the system or incorporating an image preprocessing pipeline to standardize DICOM inputs and ensure compatibility images taken from different X-ray machines. Development Plan should include risks related to false negatives (missed abnormalities) and false positives (unnecessary follow-ups). It should also include mitigation strategies such as implementing secondary review workflows.

5. Limitations of our solution includes that the AI model was trained on a specific dataset that may lack diversity in terms of patient demographics (age, gender, ethnicity) and may not perform equally well on images from underrepresented populations or different imaging machines (different brands, resolutions, or configurations). This might cause inaccuracies in model's predictions for underrepresented populations or when analyzing images from unfamiliar imaging devices. Another limitation is Module communication lacks robust retry mechanisms in case of network or hardware failures. This might cause failures in communication between system modules could result in data loss or delayed analyses. If we were given unlimited resources or had more time, we can expand the training dataset to include more diverse populations and images taken from different imaging equipments. We can also add advanced logging, monitoring, and error-recovery mechanisms to enhance fault tolerance. Our chosen design pattern is layered architecture because it is highly scalable and maintainable and easy to replace or upgrade individual layers. It is also ideal for larger systems with complex workflows since it has clear separation between system layers. For our case, our backend services has multiple AI integrations which increases complexity. Layers are interacted sequentially from top to bottom which follows the logic flow of our project. We also considered MVC as our design solution but we realized that it can lead to tight coupling between controller and view, and our project might not fit into the triangular flow of MVC.