

# Verification and Validation Report:

March 5, 2025

# 1 Revision History

Date	Version	Notes
March 01	1.0	Nanthan completed section 3
March 02	1.1	Kelly completed section 4
March 05	1.2	Patrick completed section 2, 5
March 05	1.3	Ayman completed section 7

## 2 Symbols, Abbreviations, and Acronyms

This section records information for easy reference and aims to reduce ambiguity in understanding key concepts used in the project.

### 2.1 Table of Units

Throughout this document, SI (Système International d’Unités) is employed as the unit system. In addition to basic units, several derived units are used as described below. For each unit, the symbol is given, followed by a description of the unit and the SI name.

Symbol	Unit	SI
s	Time	Second
GB	Data	Gigabyte
MB	Data	Megabyte
LOC	Quantity	Lines of Code

### 2.2 Definitions

This subsection provides a list of terms that are used in the subsequent sections and their meanings, with the purpose of reducing ambiguity and making it easier to correctly understand the requirements:

- **Artificial Intelligence (AI) Model:** A program that analyzes datasets to identify patterns and make predictions. Used extensively in medical image analysis for automating diagnostics.
- **Convolutional Neural Network (CNN):** A deep learning algorithm that processes images by assigning weights and biases, allowing it to identify patterns and features in medical images such as chest X-rays.
- **Detection Transformer (DETR):** A transformer-based neural network model for object detection. It uses an encoder-decoder transformer architecture to directly predict bounding boxes and class labels from an image, simplifying the detection process.

- **DICOM (Digital Imaging and Communications in Medicine):** The international standard for medical images, defining formats for image exchange that ensure clinical quality.
- **Containerized Application:** A portable version of an application that can be run on a container run-time, such as Docker.
- **Machine Learning (ML):** A subset of AI focusing on using data and algorithms to mimic human learning, improving accuracy over time.
- **Picture Archiving and Communication System (PACS):** A system for acquiring, storing, transmitting, and displaying medical images digitally, providing a filmless clinical environment.
- **PHI:** Personal Health Information - Private and confidential data that must be protected under the HIPPA act.
- **HIPPA:** Health Insurance Portability and Accountability Act, a set of standards protecting sensitive health information from disclosure without patient's consent.
- **AWS - Amazon Web Services:** A public cloud provider, offering all HIPAA-compliance cloud services that helps Neuralanalyzer host, manage, and scale our application.
- **AWS ECS:** AWS Elastic Cloud Service: - An AWS managed service for managing and maintaining application containers at run-time.
- **AWS ECR:** AWS Elastic Container Registry - An AWS managed service for storing and managing container images.
- **AWS Fargate:** An AWS managed service for running containerized applications.
- **AWS Cognito:** An AWS managed service for authentication logic, handling user and password management.
- **React:** A web front-end framework, written in Javascript.
- **Flask:** An HTTP-based server framework, written in Python.

- **Finite State Machine (FSM):** A computation model that simulates sequential logic using state transitions, applied in processes like user authentication and backend workflows.
- **ROC Curve (Receiver Operating Characteristic Curve):** A graph that shows the performance of a classification model by plotting the true positive rate against the false positive rate at various threshold levels.
- **Service-Level Agreement (SLA):** Defines the guaranteed uptime of the system, such as ensuring the availability of the AI service for 99.99% of operational hours.
- **Software as a Medical Device (SaMD):** Software classified as a medical device under regulatory frameworks, such as those defined by the Food and Drugs Act.
- **TorchXRAYVision:** An open-source library for classifying diseases based on chest X-ray images, offering pre-trained models to accelerate the development process.
- **X-ray:** A form of high-energy electromagnetic radiation used in medical imaging to produce images of the inside of the body, enabling the diagnosis of conditions through radiographic film or digital detectors.
- **MIT License:** An open-source software license that allows for the free use, modification, and distribution of software.
- **Training Data:** Refers to the dataset of labeled chest X-ray images used to train the AI model. In this project, the dataset size is approximately 471.12 GB.

## 2.3 Abbreviations and Acronyms

Symbol	Description
SRS	Software Requirements Specification
AI	Artificial Intelligence
CNN	Convolutional Neural Network
DICOM	Digital Imaging and Communications in Medicine
DETR	Detection Transformer
ViT	Vision Transformer
VnV	Verification and Validation
ML	Machine Learning
PACS	Picture Archiving and Communication System
SaMD	Software as a Medical Device
ROC	Receiver Operating Characteristic Curve
SLA	Service-Level Agreement
FR	Functional Requirement
NFR	Non-Functional Requirement
FSM	Finite State Machine
CXR	Chest X-Ray Project
POC	Proof of Concept
TM	Theoretical Model
AWS	Amazon Web Services
ECS	Elastic Container Service
ECR	Elastic Container Registry
CI/CD	Continuous integration and Continuous deployment
HTTP	Hypertext Transfer Protocol

# Contents

<b>1</b>	<b>Revision History</b>	<b>i</b>
<b>2</b>	<b>Symbols, Abbreviations, and Acronyms</b>	<b>ii</b>
2.1	Table of Units . . . . .	ii
2.2	Definitions . . . . .	ii
2.3	Abbreviations and Acronyms . . . . .	v
<b>3</b>	<b>Functional Requirements Evaluation</b>	<b>1</b>
3.1	FR1 . . . . .	1
3.2	FR2 . . . . .	1
3.3	FR3 . . . . .	1
3.4	FR4 . . . . .	2
3.5	FR5 . . . . .	2
3.6	FR6 . . . . .	2
3.7	FR7 . . . . .	3
3.8	FR8 . . . . .	3
3.9	FR9 . . . . .	4
3.10	FR10 . . . . .	4
3.11	FR11 . . . . .	4
3.12	FR12 . . . . .	4
3.13	FR13 . . . . .	5
<b>4</b>	<b>Nonfunctional Requirements Evaluation</b>	<b>5</b>
4.1	Look and Feel . . . . .	5
4.2	Usability and Humanity . . . . .	6
4.3	Performance . . . . .	7
4.4	Operational and Environmental . . . . .	8
4.5	Security and Privacy . . . . .	9
4.6	Maintainability and Support . . . . .	10
4.7	Cultural . . . . .	10
4.8	Legal . . . . .	11
4.9	Health and Safety . . . . .	12
<b>5</b>	<b>Comparison to Existing Implementation</b>	<b>12</b>
5.1	Existing Projects . . . . .	12
5.2	Project Structure . . . . .	13

5.3	Datasets . . . . .	14
5.4	Architectures . . . . .	15
5.5	Summary . . . . .	16
<b>6</b>	<b>Unit Testing</b>	<b>17</b>
<b>7</b>	<b>Changes Due to Testing</b>	<b>17</b>
<b>8</b>	<b>Automated Testing</b>	<b>18</b>
<b>9</b>	<b>Trace to Requirements</b>	<b>18</b>
<b>10</b>	<b>Trace to Modules</b>	<b>18</b>
<b>11</b>	<b>Code Coverage Metrics</b>	<b>18</b>

## List of Tables

5	Trace to Requirements Table . . . . .	18
6	Trace to Modules Table . . . . .	18

## List of Figures



## 3 Functional Requirements Evaluation

### 3.1 FR1

- Description: The system shall accept chest X-ray images as input from authorized users.
- Type: Manual
- Verification: Calling HTTP Server Module API to get a presigned URL for uploading the image with and without authorization token.
- Validation: Presigned URL doesn't return if the user is not authorized.
- Result: Pass

### 3.2 FR2

- Description: The system shall enable users to input additional patient symptoms, such as cough, chest pain, or fever
- Type: Automated
- Verification: Execute integration tests on the Medical Record Management Module to include additional patient symptoms.
- Validation: Check if the additional patient symptoms are included in the medical record from the Data Persistent Module.
- Result: Pass

### 3.3 FR3

- Description: The system shall analyze chest X-ray images to detect the presence or absence of specific diseases with an accuracy of 85% or higher.
- Type: Automated
- Verification: Execute Validation Script on the Disease Detection Module to check the accuracy of the AI model.

- Validation: Check if the accuracy of the AI model is 85% or higher.
- Result: Pass

### **3.4 FR4**

- Description: The system shall indicate whether a patient's condition has improved, worsened, or remained stable between scans.
- Type: Automated
- Verification: Execute integration tests on the Disease Progression Module, with pre-defined patient conditions.
- Validation: Validated the output of the Module to include improved, worsened, or remained stable conditions.
- Result: Pass

### **3.5 FR5**

- Description: The system shall generate visual aids by highlighting affected areas on the chest X-ray images.
- Type: Manual
- Verification: Manually create a medical record on the X-Ray Report View Module.
- Validation: Validated the output of the Module to include highlighted affected areas.
- Result: Pass

### **3.6 FR6**

- Description: The system shall produce a structured, human-readable report summarizing key findings, disease detection results, and progression status.
- Type: Automated

- Verification: Execute integration tests on the Internal Report Generation Service, with pre-defined patient conditions.
- Validation: Validated the output of the Service to include key findings, disease detection results, and progression status.
- Result: Pass

### **3.7 FR7**

- Description: The system shall store patient data, including images and reports, in a secure database for future reference.
- Type: Automated
- Verification: Execute integration tests on Medical Record Management Module to create a medical record with images and reports.
- Validation: Validated that all medical records, findings, and X-Ray images are stored in the Data Persistent Module.
- Result: Pass

### **3.8 FR8**

- Description: The system shall provide alerts to clinicians if significant changes in a patient's condition are detected between scans.
- Type: Manual
- Verification: Manually create a medical record on the Disease Progression Module with significant changes in a patient's condition.
- Validation: Validated that the drastic changes in a patient's condition are presence on the UI of Patient List View Module and Patient Overview Module.
- Result: Pass

### **3.9 FR9**

- Description: The system shall allow healthcare professionals to adjust treatment plans based on the X-ray analysis results
- Type: Manual
- Verification: Manually edit a medical record on the X-Ray Report View Module with adjusted treatment plans.
- Validation: Validated that the adjusted treatment plans are presence on the UI of Patient List View Module and Data Persistent Module.
- Result: Pass

### **3.10 FR10**

- Description: The system shall display confidence levels for disease detection and progression analysis results through an intuitive user interface that requires minimal training to operate.
- Type: Manual
- Verification: Manually create a medical record on the X-Ray Report View Module with known disease confidence levels.
- Validation: Validated that the confidence levels are accurately shown on the UI of X-Ray Report View Module.
- Result: Pass

### **3.11 FR11**

- CANCELLED

### **3.12 FR12**

- Description: The system shall create a new copy of a patient's X-ray before running the AI model for analysis.
- Type: Automated

- Verification: Execute integration tests on the Medical Record Management Module to make a new record.
- Validation: Validated that patient's X-ray is already copied into the Data Persistent Module before running through the Disease Detection and Progression Module.
- Result: Pass

### 3.13 FR13

- Description: The system shall support regular updates to the AI model to incorporate new data and improve accuracy over time.
- Type: Automated
- Verification: Execute integration tests on the AI Model Update Module to ensure updates are applied correctly.
- Validation: Validate that the AI model's accuracy improves over time with new data.
- Result: Pass

## 4 Nonfunctional Requirements Evaluation

### 4.1 Look and Feel

#### 1. NFR-LF1

Initial State: The system must be installed and accessible.

Input/Condition: Adjust window/level settings. Navigate between other tabs. Apply basic tools (zoom, pan, measure). Export or save images from the viewer.

Expected Output/Result: The system responds to user interactions.

Actual Output: The system responds to user interactions promptly and correctly.

Result: Pass

## 2. NFR-LF2

Initial State: The system under test is installed and accessible on the target device. The monitor brightness is set to a comfortable level to ensure consistency during the test.

Input/Condition: User interaction with the interface for typical workflows.

Expected Output/Result: Color contrast ratio is 4.5:1 or higher for text. Font sizes at least 12-14 points for normal text.

Actual Output: All system interface showed valid color contrast ratio and appropriate font size

Result: Pass

## 4.2 Usability and Humanity

### 1. NFR-UH1

Initial State: User accounts and access credentials are prepared for all healthcare professionals participating in the test. Any required files or images are pre-loaded and accessible for the tasks.

Input/Condition: Test user given a list of specified tasks and user accounts.

Expected Output/Result: User successfully login each user account and perform: upload a chest X-ray image to the system, view analysis results or processed reports, adjust image settings (e.g., window/level, zoom), and export the analysis or report to a local folder.

Actual Output: User can login as doctor / patient and upload chest X-ray image to the system, view analysis results or processed reports.

Result: Fail. Not able to adjust image settings (e.g., window/level, zoom), and export the analysis or report to a local folder.

### 2. NFR-UH2

Initial State: List of user accounts and access credentials are prepared. Full access to the user interface. Interface elements (buttons, icons, menus) are accessible for interaction.

Input/Condition: Automated script with a list of interactions specified to interact with the tabs and other listed functionalities.

Expected Output/Result: Log whether tool-tips and help content were displayed correctly for each element tested.

Actual Output: No tool bar / help content was displayed

Result: Fail

### 4.3 Performance

#### NFR-PR1

Initial State: User is logged in. Standard chest X-ray images are available for testing in PNG/JPG format (depending on system requirements).

Input/Condition: Test user interacts with the system to upload an X-Ray image.

Expected Output/Result: AI-analyzed results of the X-Ray image is displayed in user interface within 1 minute of uploading the X-ray image.

Actual Output: Disease probabilities was produced by AI model in user interface along with summaries.

Result: Pass

#### NFR-PR2

Initial State: A monitoring system is set up to log availability metrics, such as downtime events, mean time to recovery, and uptime percentage.

Input/Condition: Server outages, network issues, high number of concurrent users.

Expected Output/Result: Uptime percentage; downtime events (log the number, duration, and cause of any disruptions or downtime); mean time to recovery.

Actual: Uptime percentage of the system is above 99

Result: Pass

#### NFR-PR3

Initial State: All dependencies (database, network, storage, and image processing engine) are configured and operational. Monitoring tools are available to track system performance, including CPU and memory usage, image processing times.

Input/Condition: Collections of 20 identical or varied images, multiple user accounts. Baseline for each image processing is 20 seconds.

Expected Output/Result: Processing time (float), system resource utilization (float), images uploading success/fail (boolean).

Actual Output: Image uploading function is working. Processing time is within 20 seconds for all images. System resource utilization showed 15

Result: Pass.

## 4.4 Operational and Environmental

### 1. NFR-OE1

Initial State: The AI system and the PACS are both connected to the same hospital network and the DICOM configuration for both systems is correctly set.

Input/Condition: AI-processed results including annotations or diagnostic report as an image (in DICOM format). Metadata with patient ID (string) and instance number (string).

Expected Output/Result: Boolean indicating whether the PACS successfully stores the AI-processed results with correct format (annotated image or report). The stored result is associated with the correct patient ID and instance number in the PACS.

Actual Output: The AI system was not able to connect to a hospital network.

Result: Fail

### 2. NFR-OE2

Initial State: The network is configured normally with minimal latency and no packet loss initially. Necessary patient data and chest X-ray studies are available for retrieval.



Input/Condition: Retrieve and process a chest X-ray image, store the processed result, and send the data to external.

Expected Output/Result: Latency is within a reasonable time (190ms 220ms), logs include latency and packet loss statistics.

Actual Output: Logs showing that latency is around 190ms and no packet loss detected.

Result: Pass

## 4.5 Security and Privacy

### NFR-SR1

Initial State: AES-256 encryption libraries (such as OpenSSL or Cryptography in Python) are configured for the system. Patient data (including DICOM images and reports) is stored on the system or transmitted over the network.

Input/Condition: X-ray image and diagnostic report with patient ID and instance number. Secret key and initialization vector (IV) for AES-256 encryption.

Expected Output/Result: All patient data, including images and reports, is encrypted using AES-256 both during storage and transmission.

Actual Output: Patient data was not encrypted with AES-256 standard.

Result: Fail

### NFR-SR2

Initial State: The AI system has different user accounts with different assigned roles. Each role has specific permissions defined.

Input/Condition: Role 1: Doctor view diagnostic results and patient images. Role 2: Administrator—can manage user accounts and system configurations (user\_radiologist, user\_admin).

Expected Output/Result: Doctor should be able to view diagnostic reports and images and should not be able to manage accounts. Administrator should be able to create and manage user accounts and should not have access to patient diagnostic reports.

Actual Output: Administrator role is not defined. Doctor is able to view diagnostic reports and not be able to manage accounts.

Result: Fail.

## 4.6 Maintainability and Support

### 1. NFR-MS1

Initial State: The AI system code-base is deployed in a version-controlled environment. Documentation for each module (e.g., README files, API references, inline comments) is available in the repository.

Input/Condition: All modules within the code repository are accessible.

Expected Output/Result: All key functionalities (e.g., data ingestion, inference, reporting) are encapsulated in separate, well-defined modules. Modules can function independently with minimal coupling. Dependencies between modules are well-documented.

Actual Output: All modules are separated and follows SOLID principles. Modules are self-contained that can be modified without breaking other parts.

Result: Pass

### 2. NFR-MS2

Initial State: Code repository contains the entire AI system, automated testing framework (e.g., JUnit, Pytest) is installed and configured in the project. Tests include unit tests, integration tests, system tests, and end-to-end tests. Code coverage tool (e.g., coverage.py, JaCoCo) is integrated.

Input/Condition: Unit tests, integration tests, and end-to-end tests.

Expected Output/Result: Code coverage log.

Actual Output: Code coverage log included as part of CI/CD pipeline which will be triggered everytime there's a new commit.

Result: Pass

## 4.7 Cultural

### NFR-CR1

Initial State: The language options (English and French) are available in the settings menu.

Input/Condition: Switch language between English and French. Generate report in English and French.

Expected Output/Result: No untranslated or misaligned content should appear.

Actual Output: System does not show French content.

Result: Fail

## 4.8 Legal

### 1. NFR-LR1

Initial State: HIPAA and PIPEDA documentation are in place, development artifacts and source code are available.

Input/Condition: Development artifacts, system design documentation.

Expected Output/Result: System design aligns with the requirements of HIPAA (US) and PIPEDA (Canada). If not, all software risks are identified, evaluated, and mitigated.

Actual output: System design fully followed HIPAA (US) and PIPEDA (Canada).

Result: Pass

### 2. NFR-LR2

Initial State: ISO 13485 documentation is in place, development artifacts are available (verification and validation plans, risk analysis reports, requirements, design documents, test plans, etc.).

Input/Condition: Development artifacts, processes (software lifecycle management process, design review process, etc.).

Expected Output/Result: All development processes align with the requirements of ISO 13485 and if not, all software risks are identified, evaluated, and mitigated.

Actual Output: Development artifacts and processes follow ISO 13485 guidelines.

Result: Pass

## 4.9 Health and Safety

### 1. NFR-HS1

Initial State: User roles are configured within the system, including Radiologist, Clinician, etc.

Input/Condition: AI-generated report, Radiologist user account.

Expected Output/Result: Logs including radiologist’s review action (confirm/reject) with a timestamp.

Actual Output: Logs including doctor’s action were shown in docker as the system is running.

Result: Pass

### 2. NFR-HS2

Initial State: The AI system is deployed and accessible to radiologists, clinicians, and other users.

Input/Condition: Access the AI-generated diagnostic report and view the user interface displaying AI results.

Expected Output/Result: Disclaimer is prominently displayed and easy to read.

Actual Output: No disclaimer is displayed.

Result: Fail

## 5 Comparison to Existing Implementation

### 5.1 Existing Projects

- [harrisonchiu/xray](#)
- [N8THEPL8/ChestLenseAI](#)
- [PLAN-Lab/CheXRelFormer](#)

## 5.2 Project Structure

Project Name	FrontEnd	BackEnd	Data Base	Cloud	Deployment
<b>Our Capstone</b>	React.js	Flask API	Amazon S3	AWS	Docker
xray	N/A	N/A	N/A	N/A	Jupyter Notebook
ChestLenseAI	HTML & CSS	Flask API	N/A	N/A	Python
CheXRelFormer	N/A	N/A	N/A	N/A	Shell Command

### 5.3 Datasets

Project Name	Dataset	Size	Classes	Link
<b>Our Capstone</b>	MIMIC-CXR-JPG 2.0.0	557.6 GB	9: Lung Opacity, Pleural Effusion, Atelectasis, Enlarged Cardiac Silhouette, Pulmonary Edema/Hazy Opacity, Pneumothorax, Consolidation, Fluid Overload/Heart Failure, Pneumonia. <b>Progress:</b> No Change, Improved, Worsened.	<a href="https://physionet.org/content/mimic-cxr-jpg/2.0.0/">https://physionet.org/content/mimic-cxr-jpg/2.0.0/</a>
<b>xray</b>	Chest-xray14	42.0 GB	14: Atelectasis, Cardiomegaly, Consolidation, Edema, Effusion, Emphysema, Fibrosis, Hernia, Infiltration, Mass, Nodule, Pleural Thickening, Pneumonia, Pneumothorax.	<a href="https://nihcc.app.box.com/v/ChestXray-NIHCC/folder/37178474737">https://nihcc.app.box.com/v/ChestXray-NIHCC/folder/37178474737</a>
<b>ChestLenseAI</b>	MIMIC-CXR-JPG 2.0.0	557.6 GB	6: Atelectasis, Cardiomegaly, Consolidation, Edema, No Finding, Pleural Effusion.	<a href="https://physionet.org/content/mimic-cxr-jpg/2.0.0/">https://physionet.org/content/mimic-cxr-jpg/2.0.0/</a>
<b>CheXRelFormer</b>	MIMIC-CXR-JPG 2.0.0, MIMIC-III 1.4	557.6 GB, 6.2 GB	3: No Change, Improved, Worsened	<a href="https://physionet.org/content/mimic-cxr-jpg/2.0.0/">https://physionet.org/content/mimic-cxr-jpg/2.0.0/</a>

## 5.4 Architectures

Project Name	Neural Network	Configuration	Link (graph) (paper)
Our Capstone	DETR	MLP	<a href="https://viso.ai/wp-content/uploads/2024/02/DETR-Architecture.jpg">https://viso.ai/wp-content/uploads/2024/02/DETR-Architecture.jpg</a> <a href="https://arxiv.org/pdf/2005.12872">https://arxiv.org/pdf/2005.12872</a>
xray	ResNet	ResNet-50	<a href="https://i.ytimg.com/vi/woEs7UCaITo/maxresdefault.jpg">https://i.ytimg.com/vi/woEs7UCaITo/maxresdefault.jpg</a> <a href="https://arxiv.org/pdf/1512.03385">https://arxiv.org/pdf/1512.03385</a>
ChestLenseAI	DenseNet	DenseNet-121	<a href="https://pytorch.org/assets/images/densenet1.png">https://pytorch.org/assets/images/densenet1.png</a> <a href="https://arxiv.org/pdf/1608.06993">https://arxiv.org/pdf/1608.06993</a>
CheXRelFormer	ViT	MLP	<a href="https://www.researchgate.net/publication/383905431/figure/fig3/AS:11431281290331182@1731595235002/Structure-of-the-backbone-PVTv2.ppm">https://www.researchgate.net/publication/383905431/figure/fig3/AS:11431281290331182@1731595235002/Structure-of-the-backbone-PVTv2.ppm</a> <a href="https://arxiv.org/pdf/2106.13797">https://arxiv.org/pdf/2106.13797</a>

## 5.5 Summary

Project Name	Summary	Reference Project	Paper
<b>Our Capstone</b>	This project encompasses the entire product lifecycle, integrating Frontend, Backend, Cloud Infrastructure, and a CI/CD pipeline. It is designed to support the detection of diseases, track the progression of conditions over time, and generate AI-driven diagnostic reports for healthcare professionals.	<a href="https://github.com/McMasterAIHLab/CheXDetector">https://github.com/McMasterAIHLab/CheXDetector</a>	<a href="https://papers.miccai.org/miccai-2024/paper/3269_paper.pdf">https://papers.miccai.org/miccai-2024/paper/3269_paper.pdf</a>
xray	This project only contains backend structure for the AI model, which uses a smaller dataset compared to other projects and uses the resnet-50 network for disease classification.	<a href="https://github.com/LalehSeyyed/CheXclusion">https://github.com/LalehSeyyed/CheXclusion</a>	<a href="https://arxiv.org/pdf/2003.00827v2">https://arxiv.org/pdf/2003.00827v2</a>
ChestLenseAI	This project uses a minimal Frontend and applies the DenseNet-121 pre-trained model to detect diseases in chest X-ray images.	<a href="https://github.com/LaurentVeyssier/Chest-X-Ray-Medical-Diagnosis-with-Deep-Learning/tree/main">https://github.com/LaurentVeyssier/Chest-X-Ray-Medical-Diagnosis-with-Deep-Learning/tree/main</a>	<a href="https://arxiv.org/pdf/1711.05225">https://arxiv.org/pdf/1711.05225</a>
CheXRelFormer	This project is purely focused on the backend ViT model, it does not use a pre-trained network for disease progress between two x-ray images. It is a very advanced (PhD) project which is developed with no prior related research.	N/A	N/A



## 6 Unit Testing

## 7 Changes Due to Testing

To run automated tests within a developer's local environment, a developer can execute a command to build the project. This process ensures that all dependencies are correctly set up and runs the automated tests for both Python and React.js.

If a developer wants to run only the tests, they can navigate to the appropriate directory where the tests are located and execute the relevant test command. For Python, this typically involves running `pytest` or `Flake8`, while for React.js, `Jest` is used.

We have tests that are executed using AWS services, with results and logs stored in Amazon S3 for easy access and review. This setup ensures that test results are centralized and accessible across for the team. The automation helps maintain code quality, catch errors early, and streamline the deployment process.

From the repository's perspective, tests are executed using Github CI/CD to maintain code quality and stability. Both linting and compilation are performed using the same commands that a developer would execute in their local environment. Specifically, we use `Flake8` for Python linting and `Jest` for testing React.js components.

These tests are triggered when a new pull request is made to the main branch. If any test fails, merging is blocked until the issues are resolved. This ensures that only verified, high-quality code is integrated. Additionally, a compiler workflow runs after merging into the main branch to detect any errors or unintended changes in code behavior.

If any test or workflow fails, the detailed logs provide valuable insights into the reason for failure. This allows developers to quickly diagnose and address issues, preventing regressions and ensuring a smooth development process. By enforcing these automated checks, we maintain code consistency, reduce bugs, and improve overall project reliability.

## 8 Automated Testing

## 9 Trace to Requirements

For additional information about our functional requirements, please refer to our [SRS](#).

Column 1	Column 2
Cell 1	Cell 2
Cell 3	Cell 4

Table 5: Trace to Requirements Table

## 10 Trace to Modules

For additional information about our modules, please refer to our [MG](#).

Column 1	Column 2
Cell 1	Cell 2
Cell 3	Cell 4

Table 6: Trace to Modules Table

## 11 Code Coverage Metrics

## References

## Appendix — Reflection

The information in this section will be used to evaluate the team members on the graduate attribute of Reflection.

1. What went well while writing this deliverable?
2. What pain points did you experience during this deliverable, and how did you resolve them?
3. Which parts of this document stemmed from speaking to your client(s) or a proxy (e.g. your peers)? Which ones were not, and why?
4. In what ways was the Verification and Validation (VnV) Plan different from the activities that were actually conducted for VnV? If there were differences, what changes required the modification in the plan? Why did these changes occur? Would you be able to anticipate these changes in future projects? If there weren't any differences, how was your team able to clearly predict a feasible amount of effort and the right tasks needed to build the evidence that demonstrates the required quality? (It is expected that most teams will have had to deviate from their original VnV Plan.)